


GEMINI

Uroflowmeters designed with users in mind

Affordability meets functionality



Better uroflowmeters for your practice



Self Calibrates

Automatically calibrates when a beaker is placed, ensuring accurate measurements with various beaker types



Easily Adjustable

One-hand adjustable stand easily accommodates standing patient or placed under commode



Automatic

Recording starts automatically when your patient starts voiding



Waterproof

Waterproof weight cell improves hygiene and simplifies cleaning

Danflow Wave CU

Uroflowmetry that's quick, easy, accurate, and provides real-time reporting without the need for a laptop.

- Thermal printer provides immediate zebra printouts of the test for analysis
- Ideal for fast, standalone uroflowmetry testing



Danflow Wave

Smart uroflowmeter that wirelessly connects to your laptop and offers advanced options.

- Customizable recording options and comprehensive reports with connection to patient information that display graph and nomogram
- Ideal for advanced and efficient uroflowmetry testing



Danflow Wave with Laptop and Integrated Bladder Scanner

Includes all the features of the Danflow Wave, with post-void residual (PVR) added to provide convenience for practices.

- Integrated bladder scanner and laptop enable uroflowmetry and PVR measurement at the same time
- Ideal for streamlined workflows and reporting with uroflowmetry and PVR consolidated



Choose the uroflowmeter that best meets your practices' needs

FEATURES		DANFLOW WAVE CU	DANFLOW WAVE	DANFLOW WAVE WITH LAPTOP & INTEGRATED BLADDER SCANNER
Automatic calibration & recording		Yes		
Adjustable stand		Yes, one-handed		
Waterproof weight cell		Yes		
Report Printing		Via built-in thermal printer Zebra printout	Via wirelessly connected laptop (not included) Comprehensive	Via wirelessly connected laptop (included) Comprehensive
EMR		-	Optional HL7, DICOM, save to networked folder	
Warranty	Components	3 years*		
	Computer	-	-	3 years
	Accessories	3 years*	-	3 years*
	Downtime	Advanced Exchange		
	Remote Service/Support	Limited	3 years	

*Probe has 2-year warranty and Accidental Drop Coverage (ACD) not available

Get more info at:
geminimedtech.com/uroflowmeters



Gemini Medical Technologies LLC

Phone: +1.833.943.6464
 Email: info@geminimedtech.com
geminimedtech.com

123 NW 13th Street, Suite 212 Boca Raton, FL 33432